





FAQ - JENOPTIK GRYPHAX® which 3rd party devices from PRIOR Scientific are supported?

The Prior device driver for GRYPHAX application uses PRIOR SDK / DLL library for Windows OS.

In general, GRYPHAX will handle and control the devices, which are available and made visible by the PRIOR scientific SDK version 1.86.

The XYZ-Stage tool of GRYPHAX is compatible to the following hardware from PRIOR Scientific:


PRIOR ProScan I, II, III (H29/H30/H31) and OptiScan (ES9/ES10/ES11) controllers as well as PRIOR motorized scanning Stages (H1XXXX) and focus motors for Z-axis (motor stepper).

Please note, we test the PRIOR ProScan III (H31) controller together with X,Y motorized scanning stage (H101A) and focus motor Z-axis (motor stepper) only.

(Only certain tests with piezo Z-drive as well).

Other PRIOR equipment are not tested so far.


Detailed information about tested PRIOR controller:







GRYPHAX XY Scan tool during automatic image stitching supported by PRIOR motorized stage:



Download latest application

Free download of <u>GRYPHAX</u> after registration: <u>https://www.ienoptik.com/progres-gryphax-software</u>

Latest GRYPHAX software manual within all new feature and tools is download able from the link below: https://www.jenoptik.com/-/media/websitedocuments/optics/progres/software-gryphax/software-manual.pdf

Learn more about JENOPTIK GRYPHAX® software and the easy workflows and tools.

Watch our video tutorials.

HELPFUL? We appreciate your feedback.



We looking forward to hearing from you and remain sincerely,

Your JENOPTIK GRYPHAX® Team

